FRONT

care
 GROUP
 EDO CHAIN

PATIENT INFORMATION LEAFLET

Device Type: Intraocular Lens (Phakic)

Introduction

This patient information leaflet is designed to help you decide whether you want to have an Implantable Phakic Contact Lens (IPLC V2.0) placed in one or both eyes to correct your nearsightedness (myopia) or farsightedness (hyperopia). It is important to understand the benefits and possible risks before making your decision. Please read this leaflet carefully and discuss any questions with a doctor trained in IPLC V2.0 surgery.

What is Nearsightedness and Farsightedness?

NORMAL EYE

Light focuses on the retina.
 Vision is clear.

Normal Eye

NORMAL EYE:

Your eye works like a camera. You see objects because your eye focuses light into images. Both the clear front layer of your eye, called the cornea, and your natural crystalline lens focuses rays of light onto the back surface of the eye, called the retina.

NEARSIGHTED EYE

Light focuses in front of the retina. Vision is blurry.

Nearsighted Eye

NEARSIGHTED EYE:

Nearsightedness is a focusing error that causes blurrer vision at distance than at near. Light from an object in the distance focuses in front of the back surface of your eye (retina) responsible for capturing images, rather than directly on the retina.

FARSIGHTED EYE

Light focuses behind the retina. Vision is blurry.

Farsighted Eye

FARSIGHTED EYE:

Farsightedness (hyperopia) is a common vision condition in which you can see distant objects clearly, but objects nearby may be blurry.

Glasses and contact lenses can help your eye focus light properly on the retina. The IPLC V2.0 can also help your eye focus light properly on the retina.

How Does IPLC V2.0 Correct Nearsightedness and farsightedness?

CORRECTION OF VISION AFTER IPLC V2.0 SURGERY

Eye After Treatment

A diopter (D) is a unit to measure how well your eye can focus. It shows the level of nearsightedness, farsightedness, or astigmatism you have. IPLC V2.0 surgery can treat nearsightedness from +0.5 D to -15.0 D and farsightedness from +0.5 D to +15.0 D. If your vision falls within these ranges, this surgery may help you see clearly without glasses or contact lenses.

The IPLC V2.0 monofocal lens does not correct astigmatism. Astigmatism is a common vision problem that causes blurred or uneven vision. It happens when the front surface of the eye (cornea) is not evenly shaped. Because of this, light does not focus properly inside the eye, and some parts of what you see may look clearer than others.

NEARSIGHTED AND ASTIGMATIC EYE

Nearsighted and Astigmatic Eye

IPLC V2.0 Toric / IPLC V2.0 Presbyopic: Toric can correct any astigmatism you have up to +1.0 D to +10D. You are requested to consult your doctors to choose appropriate suitable options for your eyes.

Choosing the best lens for your vision

There are different types of lenses available. Your eye doctor will help you choose the one that is best for you. Other treatment options include glasses with thick lenses or contact lenses instead of a permanent lens implant. Talk to your eye doctor about the benefits and possible risks of each option.

Other Treatments to Correct Nearsightedness and Farsightedness ?

Other treatments to correct nearsightedness and farsightedness include eyeglasses, contact lenses, and refractive surgical procedures. Surgical options include Photorefractive Keratectomy (PRK), Laser-Assisted In Situ Keratomileusis (LASIK), refractive lens exchange (RLE), and phakic intraocular lens implantation (e.g., IPLC V2.0). However, these procedures may not meet the visual requirements for certain professions, such as military service.

These procedures either reshape the cornea or place a lens inside the eye. PRK and LASIK use a laser to reshape the shape of the cornea so the eye can focus better. In LASIK, a thin flap of the cornea is lifted, the laser reshapes the tissue underneath, and the flap is replaced in place. Phakic intraocular lens implantation, like IPLC V2.0, adds an artificial lens inside the eye without removing your natural lens.

Benefit And Risks of the IPLC V2.0

IPLC V2.0 surgery can safely correct nearsightedness between -0.50 to -15.0D and farsightedness between +0.50 to +15.0 D. If your eyeglasses prescription is in these ranges, the IPLC V2.0 may make your vision clearer without eyeglasses or contact lenses correction better.

Implanting IPLC V2.0 requires eye surgery, which carries potential risks. Please read this leaflet carefully and talk to your doctor about the possible risks.

■ Few of the risks are as below

■ Additional surgery
 Sometimes, a second surgery may be needed to remove, replace, or adjust the IPLC V2.0. It may need to be removed during cataracts surgery. Additional surgery might also be needed to reduce eye pressure.

■ Cataract formation
 After Phakic IOL surgery, there is an increased risk that the natural lens of your eye may become cloudy, leading to cataracts that could require surgery. The longer the IPLC V2.0 stays in the eye, the higher this risk becomes. Regular eye exams are important to check for cataracts. Older patients and those with stronger nearsightedness or farsightedness may have a higher risk. The long-term risk beyond 7 years is not known. If cataract surgery is needed, the IPLC V2.0 and the cloudy lens are removed, and a new lens is implanted, likely the modified cataract surgery.

■ Loss of Best Corrected Vision (BCVA)
 The most common reason was clouding of the natural lens (cataract). Vision may get better in these eyes after cataract surgery.

■ Raised pressure inside the eye and damage to the optical nerve caused by increased pressure in eye (Glaucoma)
 Normal eye pressure (intraocular pressure or IOP) can vary but is often considered to be from 10-21 millimeters of mercury (mmHg). An IOP higher than normal is called ocular hypertension or if left untreated, can cause damage to the optic nerve (glaucoma) and cause permanent vision loss.

■ Endothelial Cell Loss

Removal of the clear front part of your eye has a thin layer of special cells called endothelial cells. These cells keep the cornea clear by removing the excess fluid. Normally, the number of these cells slowly decreases with age. Eye surgery, including IPLC V2.0, can cause additional loss of these cells. If too many cells are lost, the cornea can become cloudy and affect your vision. Before surgery, your doctor will examine your eyes to check if you are suitable for IPLC V2.0. After the surgery, some ongoing loss of these cells may occur over time. If too many cells are lost, fluid can build up in the cornea (called corneal edema), which may require a cornea transplant.

■ Other Complications

Other risks associated with IPLC V2.0 surgery may include:
 ■ movement of the contact portion of the lens (iris) through a surgical wound to a position outside the eye (iris prolapse),
 ■ bleeding in the area or the natural response for reading vision (macular hemorrhage),
 ■ bleeding under the retina (subretinal hemorrhage),
 ■ increase in focusing error (astigmatism),
 ■ lifting or pulling of the retina from its natural position (retinal detachment), and
 ■ unequal pupil size (anisocoria).

Adverse reactions and complications due to, or following surgery and implantation of any implantable Phakic Contact Lens may include, but are not limited to:

- Increased astigmatism
- Refractive error
- Intraocular infection
- Vitreous loss
- Visual loss
- Anterior subcapsular cataract
- Severe intraocular inflammation
- Glaucoma
- Secondary glaucoma
- Uveitis
- Vitritis
- Hyphema
- Hyphema
- Hyphema</

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