FRONT

care
GROUP

PATIENT INFORMATION LEAFLET

Device Type: Intraocular Lens (Phakic)

Introduction


This patient information leaflet is designed to help you decide whether you want to have an Implantable Phakic Contact Lens (EYEPCL V2.0) placed in one or both eyes to correct your nearsightedness (myopia) or farsightedness (hyperopia). It is important to understand the benefits and possible risks before making your decision. Please read this leaflet carefully and discuss any questions with a doctor trained in EYEPCL V2.0 surgery.

What is Nearsightedness and Farsightedness?

NORMAL EYE

Light focuses on the retina.

Vision is clear.


NORMAL EYE:

Your eye works like a camera. You see objects because your eye focuses light into images. Both the clear front layer of your eye, called the cornea, and your natural crystalline lens focuses rays of light onto the back surface of the eye, called the retina.

NEARSIGHTED EYE:

Nearsightedness is a focusing error that causes blurrer vision at distance than at near. Light from an object in the distance focuses in front of the back surface of your eye (retina) responsible for capturing images, rather than directly on the retina.


NEARSIGHTED EYE:

Light focuses in front of the retina. Vision is blurry.

FARSIGHTED EYE:

Light focuses behind the retina. Vision is blurry.

FARSIGHTED EYE:


Farsightedness (hyperopia) is a common vision condition in which you can see distant objects clearly, but objects nearby may be blurry.

Glases and contact lenses can help your eye focus light properly on the retina. The EYEPCL V2.0 can also help your eye focus light properly on the retina.


How Does EYEPCL V2.0 Correct Nearsightedness and Farsightedness?

CORRECTION OF VISION AFTER EYEPCL V2.0 SURGERY

NEARSIGHTED EYE

Eye After Treatment



A diopter (D) is a unit to measure how well your eye can focus. It shows the level of nearsightedness, farsightedness, or astigmatism you have. EYEPCL V2.0 surgery can correct nearsightedness from -0.5 D to -15.0 D and farsightedness from +0.5 D to +15.0 D. If your vision falls within these ranges, the surgery may help you see clearly without glases or contact lenses.

The EYEPCL V2.0 monofocal lens does not correct astigmatism. Astigmatism is a common vision problem that causes blurred or uneven vision. It happens when the front surface of the eye (cornea) is not evenly shaped. Because of this, light does not focus properly inside the eye, and some parts of what you see may look clearer than others.

NEARSIGHTED AND ASTIGMATIC EYE

Nearsighted and Astigmatic Eye

EYEPCL V2.0 Toric / EYEPCL V2.0 Presbyopic Toric can correct any astigmatism you have up to +1.0 D to +10.0. You are requested to consult your doctor to choose the appropriate suitable options for your eye.

Choosing the best lens for your vision

There are different types of lenses available. Your eye doctor will help you choose the one that is best for you. Other treatment options include glases with thick lenses or contact lenses instead of a permanent lens implant. Talk to your eye doctor about the benefits and possible risks of each option.

Other Treatments to Correct Nearsightedness and Farsightedness?

Other treatments to correct nearsightedness and farsightedness include eyeglasses, contact lenses, and refractive surgical procedures. Surgical options include Photorefractive Keratectomy (PRK), Laser-Assisted In Situ Keratomileusis (LASIK), refractive lens exchange (RLE), and phakic intraocular lens implantation (e.g., EYEPCL V2.0). However, these procedures may not meet the visual requirements for certain professions, such as military service.

These procedures either reshape the cornea or place a lens inside the eye. PRK and LASIK use a laser to change the shape of the cornea so the eye can focus better. In LASIK, a thin flap of the cornea is lifted, the laser reshapes the tissue underneath, and the flap is replaced. If you choose to have phakic intraocular lens implantation, like EYEPCL V2.0, adds an artificial lens inside the eye to correct your vision.

Benefit And Risks of the EYEPCL V2.0

EYEPCL V2.0 surgery can safely correct nearsightedness between -0.50 to +15.0 and farsightedness between +0.50 to +15.0. If your eyeglass prescription is in these ranges, the EYEPCL V2.0 may allow you distance vision without eyeglasses or contact lens correction better.

Implanting EYEPCL V2.0 requires eye surgery, which carries potential risks. Please read this leaflet carefully and talk to your doctor about the possible risks.

Few of the risks are as below

- Additional surgery**
Sometimes, a second surgery may be needed to remove, replace, or adjust the EYEPCL V2.0. It may need to be removed during cataract surgery. Additional surgery might also be needed to reduce eye pressure.
- Cataract formation**
After Phakic IOL surgery, there is an increased risk that the natural lens of your eye may become cloudy, leading to cataracts that could require surgery. The longer the EYEPCL V2.0 stays in the eye, the higher this risk becomes. Regular eye exams are important to check for cataracts. The eyes and those with stronger nearsightedness or farsightedness may have a higher risk. The long-term risk beyond 7 years is not known. If cataract surgery is needed, the EYEPCL V2.0 and the cloudy lens are removed, and a new lens is implanted, until the next refraction surgery.
- Loss of Best Corrected Vision (BCVA)**
The most common reason was clouding of the natural (cataract) vision. You may get better in these eyes after cataract surgery.
- Raised pressure inside the eye and damage to the optical nerve caused by increased pressure in eye (Glaucoma)**
Normal eye pressure (intraocular pressure or IOP) can vary but is often considered to be from 10-21 millimeters of mercury (mmHg). An IOP higher than normal is called ocular hypertension or if left untreated, can cause damage to the optic nerve (glaucoma) and cause permanent vision loss.
- Endothelial Cell Loss**
The cornea (the clear front part of your eye) has a thin layer of special cells called endothelial cells. These cells keep the cornea clear by removing extra fluid. Normally, the number of these cells slowly decreases with age. Eye surgery, including EYEPCL V2.0, can cause additional loss of these cells. If too many cells are lost, the cornea can become cloudy and affect your vision. Before surgery, your doctor will examine your eyes to check if you are suitable for EYEPCL V2.0. After the surgery, some ongoing loss of these cells may occur over time. If too many cells are lost, fluid can build up in the cornea (called corneal edema), which may require a cornea transplant.
- Other Complications**
Other risks associated with EYEPCL V2.0 surgery may include:
 - movement of the colored portion of the eye (iris) through a surgical wound to a position outside the eye (iris prolapse),
 - check for extra fluid. Normally, the number of these cells slowly decreases with age. Eye surgery, including EYEPCL V2.0, can cause additional loss of these cells. If too many cells are lost, the cornea can become cloudy and affect your vision. Before surgery, your doctor will examine your eyes to check if you are suitable for EYEPCL V2.0. After the surgery, some ongoing loss of these cells may occur over time. If too many cells are lost, fluid can build up in the cornea (called corneal edema), which may require a cornea transplant.
 - bleeding under the retina (subretinal hemorrhage),
 - increase in focusing error (accommodation),
 - lifting or pulling of the retina from its natural position (retinal detachment),
 - unequal pupil size (anisocoria).

Adverse reactions and complications due to, or following surgery and implication of any implantable Phakic Contact Lens may include:

- Increased astigmatism
- Refractive error
- Intraocular infection
- Vitreous loss
- Waisting issue
- Anterior subcapsular cataract
- Severe intraocular inflammation
- Glaucoma
- Secondary glaucoma
- Uveitis
- Vitritis
- Hyphema
- Hyphema
- Hyphema
- Under/over correction
- Non-reactive pupil
- Endothelialitis
- TASS
- Retinal Detachment
- Iritis
- Ocular macular edema
- Macular edema
- Corneal decompensation
- Corneal edema
- Pupillary block
- Iris Capture
- IOP > 30mmHg over baseline and IOP > 25 mmHg
- Significant glares and/or halos after surgery
- Decentration/rotation of the EYEPCL V2.0
- Corneal endothelial cell loss
- Lens dislocation/rotation of the EYEPCL V2.0
- Loss of best spectacle corrected visual acuity

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375 x 250 mm

- Secondary surgical intervention to remove/replace the EYEPCL V2.0
- Implantation of an EYEPCL V2.0 may result in a decrease in central endothelial cell density

Contraindications and Warnings

- Implantable Phakic Contact Lens (EYEPCL V2.0) is contraindicated in the presence of any of the following circumstances and/or conditions but not limited to the following: Patients with a history of/for clinical signs of iris, uveitis, synchysis, pigment dispersion syndrome, retinal disease (other than manifestation of myopia degeneration), chronic intraocular inflammation, macular degeneration, irregular astigmatism, keratoconus or ectatic macular edema in either eye
- Patients with diabetic retinopathy in either eye
- Patients with contact glaucoma in either eye
- Patients with history of previous intraocular surgery (including refractive surgery) in the eye to be treated
- Patients with a history of uveitis in either eye
- Patients who are amblyopic or blind in the fellow eye
- Patients who do not fall in the range of pre-treatment myopia and/or astigmatism as outlined in this
- clinical prescreen inclusion criteria
- Patients with a progressive sight-threatening disease. Patients with retinal findings associated with pathological myopia are allowed.
- Patients with a diagnosis of ocular hypertension
- Patients with contact glaucoma in either eye
- Patients with Fuchs's dystrophy or other central pathology
- Patients with cataract in either eye
- Patients with an unstable refractive error in the past 12 months (myopia and/or astigmatism)
- Patients who are pregnant or nursing
- Patients with low/normal central endothelial cell density
- Patients with Anterior chamber depth <2.8mm

Are You A Good Candidate For EYEPCL V2.0 Surgery?

Your corneal topograpy a detailed eye check-up to see if EYEPCL V2.0 surgery is suitable for you. You may be considered for this surgery if you:

- are 20 to 45 years old (for Monofocal and Toric lenses) or 40 to 60 years old (for presbyopic lenses)
- have nearsightedness between -0.5 D to -15.0 D or farsightedness between +0.5 D to +1.0 D (for EYEPCL V2.0 Monofocal and EYEPCL V2.0 Toric) and astigmatism between +1.0 D and +10.0 D (for EYEPCL V2.0 Presbyopic and EYEPCL V2.0 Presbyopic Toric)
- have a regular astigmatism (for toric lenses), if applicable
- have enough space inside your eye to safely place the lens
- have a healthy endothelium (the inner layer of your eye) with enough cells (endothelial cells) on the innermost surface of the cornea, responsible for keeping the cornea clear. If your doctor determines that your endothelial cell density is below the minimum level, you will be greater risk of swelling of your cornea (corneal edema), possibly requiring removal and replacement of your corneal transplant
- have written evidence that your nearsightedness or farsightedness has been stable for at least 1 year;
- have no corneal irregularities or scars that may affect your ability to drive after the procedure, also you must
- have no known allergies to any of the medications that your doctor may discuss will be used before, during and after the procedure
- are not pregnant or breastfeeding
- understand that at between 2 to 3 weeks before EYEPCL V2.0 surgery you will need to have holes made in the extreme outer edge of the colored portion of the eye (peripheral iridectomy) using a laser
- are willing to sign a consent form before the procedure

You and your doctor will determine if you are a suitable candidate for the EYEPCL V2.0 (monofocal, toric, presbyopic and presbyopic toric) and the frequency of follow-up visits to monitor and adjust the health of your eye.

What Should You Expect During EYEPCL V2.0 Surgery?

Before the Surgery

Before surgery, your doctor will review the medical and eye history and examine both eyes to make sure you are healthy and suitable for EYEPCL V2.0. This includes the frequency of the inner layer of your cornea (endothelium). If you wear contact with low/normal central endothelial cell density, you should not wear them for at least 2 to 4 weeks before your eye check-up and surgery. This helps your doctor get accurate measurements. Not doing this may affect your results.

After your doctor has given you any medicines you may need, eye problems, past eye surgeries, medical conditions, or allergies. Ask if you should not eat, drink, before surgery. Arrange for someone to take you home, so you should not drive right after the procedure. Your doctor will tell you when it is safe to drive again.

Two to Three Weeks before Surgery

Two to three weeks before your EYEPCL V2.0 surgery, your doctor will make two holes in the extreme outer edge of the colored portion of the eye (peripheral iridectomy) using a laser. This procedure is done to prevent the possibility of angle closure glaucoma. If the fluid flows properly from the back chamber to the front chamber of the eye to prevent a buildup of pressure within the eye after the lens is placed. Your doctor will usually apply numbing drops to the eye and make tiny incisions in the colored portion of the eye (iris) with a laser beam. Usually, this procedure does not affect your ability to drive after the procedure but check with your doctor. After the peripheral iridectomy procedure, you will be prescribed eye drops for you to use. It is important that you follow all medication instructions. Your doctor will instruct you to discontinue the use of these medications before the day of surgery.

On the Day of Surgery

On the day of surgery, eye drops will be placed in your eye to enlarge (dilate) the black part of your eye (pupil). Once your pupil is fully dilated, your doctor will put numbing drops in your eye and/or inject with numbing medication into your eye and ask you to lie on your back on the treatment table in the treatment room. Your doctor may discuss different options for numbing or relaxing you before the surgery. A small cut is made in the cornea (the clear front part of the eye), and the EYEPCL V2.0 lens is placed in the eye and centered inside your eye. The procedure usually takes about 20 to 30 minutes.

After the surgery, your doctor will apply eye drops or ointment and may cover your eye with a patch or shield for protection. The procedure is painless due to the numbing medicines. Do not drive yourself home after surgery. Arrange someone for someone to take you home in advance.

The First Day after Surgery

Your doctor will check your vision the day after surgery, including measuring the pressure in your eye. You may feel sensitive to light or feel like something is in your eye. Wearing sunglasses can help. Some discomfort or mild pain is normal, and your vision may give you difficulty to see clearly. If you have any of these symptoms, you should not use eye drops or ointment for about a week to prevent infection and reduce inflammation.

IMPORTANT: Use the eye medications as directed by your eye doctor. (Your results may depend upon you following your doctor's instructions).

DO NOT touch your eyes, especially for the first 3 to 5 days after surgery. If you notice any sudden decrease in vision, you should contact your doctor immediately. Wear sunglasses if you experience any of the following symptoms:

Long Term Care: In some cases, the EYEPCL V2.0 may need to be replaced or revised. It may be replaced if another lens is better suited for your eye or you can give you better vision. It may be revised if you develop a cataract and need cataract surgery in such cases, a new lens can be implanted to help correct your vision.

EYEPCL V2.0 is removed, you will lose your vision completely if removed. Your vision may not be the same as before surgery. It is important to follow your doctor's advice and attend all follow-up visits after surgery.

Questions To Ask Your Doctor

You may want to ask the following questions to help decide if EYEPCL V2.0 surgery for nearsightedness or farsightedness is right for you:

- What other options do I have to correct my nearsightedness or farsightedness?
- Will I have to limit my activities after surgery and for how long?
- How will be the benefits of EYEPCL V2.0 surgery and for what amount of nearsightedness and farsightedness?
- How long will I have to wait before I can have surgery in my other eye?
- How long will I have to wait before I can have surgery in my other eye?
- What vision problems might I experience if I have EYEPCL V2.0 even in 1 eye?
- How often? Could my need for eyeglasses increase over time? Could I undergo a different type of eye surgery for the correction of my vision?
- How long? EYEPCL V2.0 surgery likely to affect my need to wear eyeglasses or contacts less as 1 year after?
- Will my eyes heal differently, if injured after implantation of the EYEPCL V2.0?
- Should I have EYEPCL V2.0 surgery in my other eye?
- How long will I have to wait before I can have surgery in my other eye?
- What vision problems might I experience if I have EYEPCL V2.0 even in 1 eye?

Summary of Important Information

- EYEPCL V2.0 surgery provides a permanent correction of your nearsightedness or farsightedness as long as the EYEPCL V2.0 remains in the eye. The EYEPCL V2.0 may need to be replaced, if it is removed, you will lose this benefit, and your vision may not be the same as before surgery.
- EYEPCL V2.0 surgery does not eliminate the need for reading glasses, even if you have never worn them before.
- Quantities to follow your doctor's advice and attend all follow-up visits after surgery.
- Your vision may be stable before EYEPCL V2.0 surgery. You must provide written evidence that your nearsightedness or farsightedness has changed by no more than 0.50 D each year for at least 1 year.
- If you are pregnant or breastfeeding, you should wait before having the surgery.
- Like all surgeries, EYEPCL V2.0 has risks. It is important to read and understand all the information before deciding.
- Other options to correct your vision include glasses, contact lenses, and surgeries like RK, PRK, and LASIK.
- Before considering EYEPCL V2.0 surgery you should:
- have a complete eye examination.
- talk with at least one eye care professional about EYEPCL V2.0 surgery, especially the potential benefits, risks, and complications. You should discuss the time needed for healing after surgery.

Magnetic Resonance Compatibility: The implanted IOL does not contain any metallic components and is safe with no potential for interaction with the Magnetic Resonance (MR) field. IOL is safe for Magnet Resonance Imaging (MRI) at a magnetic field strength of 1.5 and 3.0 Tesla.

Expected benefits and necessary follow-up
The IOL will usually last for your lifetime. Should you however encounter any deterioration in your vision quality, please consult your eye doctor for further investigation and treatment. Your eye doctor will discuss with you the necessary follow-up intervals based on your medical condition.

Reporting of serious incidents

Any incident related that occurs in connection to the intraocular lens should be reported to your doctor (in the first instance), to the local representative or to the CAREGROUP SIGHT SOLUTION PRIVATE LIMITED.

About this document

Your eye doctor has performed the refractive surgery and implanted the intraocular lens (IOL) in your eye. The patient implant card and leaflet have been given to you after the surgery. Please always carry the patient implant card with you to facilitate medical care.

This document is provided to the implant card and provide you latest updates about the IOL.

Mfg. Lic. No.: MFG/MD/2020/00163
Manufactured by:
CAREGROUP SIGHT SOLUTION PRIVATE LIMITED
Block No.: 31/0/CE, Village Dabhata, Taluka - Padra,
Distt. Vadodra - 391 440, Gujarat, India
Tel No.: 1800 120 3688
customercare@cargroupindia.com
www.caregroupindia.com

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