
Introduction
This patient information leaflet is designed to help you decide whether you want to have an Implantable Phakic Contact Lens (IPL V2.0) placed in one or both eyes to correct your nearsightedness (myopia) or farsightedness (hyperopia). It is important to understand the benefits and possible risks before making your decision. Please read this leaflet carefully and discuss any questions with a doctor trained in IPL V2.0 surgery.


What Is Nearsightedness and Farsightedness?

NORMAL EYE
Light focuses on the retina. Vision is clear.


NORMAL EYE
Your eye works like a camera. You see objects because your eye focuses light into images. Both the clear front layer of your eye, called the cornea, and your natural crystalline lens focuses rays of light onto the back surface of the eye, called the retina.

NEARSIGHTED EYE
Light focuses in front of the retina. Vision is blurry.

NEARSIGHTED EYE
Your eye works like a camera. You see objects because your eye focuses light into images. Both the clear front layer of your eye, called the cornea, and your natural crystalline lens focuses rays of light onto the back surface of the eye (retina) responsible for capturing images rather than directly on the retina.

FARSIGHTED EYE
Light focuses behind the retina. Vision is blurry.


FARSIGHTED EYE:
Farsightedness (hyperopia) is a common vision condition in which you can see distant objects clearly, but objects nearby may be blurry.

Glasses and contact lenses can help your eye focus light properly on the retina. The IPL V2.0 can also help your eye focus light properly on the retina.

How Does IPL V2.0 Correct Nearsightedness and Farsightedness?

CORRECTION OF VISION AFTER IPL V2.0 SURGERY


A diopter (D) is a unit used to measure how well your eye can see. It shows the level of nearsightedness, farsightedness, or astigmatism you may have. IPL V2.0 surgery can treat nearsightedness from -0.5 D to -15.0 D and farsightedness from +0.5 D to +15.0 D. If your vision falls within these ranges, this surgery may help you see clearly without glasses or contact lenses.

The IPL V2.0 monofocal lens does not correct astigmatism. Astigmatism is a common vision problem that causes blurred or uneven vision. It happens when the front surface of the eye (cornea) is not evenly shaped. Because of this, light does not focus properly inside the eye, and some parts of what you see may look clearer than others.

NEARSIGHTED AND ASTIGMATIC EYE

IPL V2.0 Toric / IPL V2.0 Presbyopic Toric can correct any astigmatism you have up to +1.0 D to +10D. You are requested to consult your doctor to choose appropriate suitable options for your eye.

Choosing the best lens for your vision

There are different types of lenses available. Your eye doctor will help you choose the one that is best for you. Other treatment options include glasses with thick lenses or contact lenses instead of a permanent lens implant. Talk to your eye doctor about the benefits and possible risks of each option.

Other Treatments to Correct Nearsightedness and Farsightedness ?
Other treatments to correct nearsightedness and farsightedness include eyeglasses, contact lenses, and refractive surgical procedures. Surgical options include Photorefractive Keratectomy (PRK), Laser-Assisted in Situ Keratomileusis (LASIK), refractive lens exchange (RLE), and phakic intraocular lens implantation (e.g., IPL V2.0). However, these procedures may not meet the visual requirements for certain professions, such as military service.

These procedures either reshape the cornea or place a lens inside the eye. PRK and LASIK use a laser to reshape the shape of the cornea so the eye can focus better. In LASIK, a thin flap of the cornea is lifted, the laser reshapes the tissue underneath, and the flap is put back in place to heal. Phakic intraocular lens implantation, like IPL V2.0, adds an artificial lens inside the eye without removing any natural lens.

Benefit And Risks of the IPL V2.0

IPL V2.0 surgery can safely correct nearsightedness between -0.5D to -15D and farsightedness between +0.5D to +15D. If your eyeglass prescription is in these ranges, the IPL V2.0 may make your distance vision without eyeglasses or contact lenses. Implanting IPL V2.0 requires eye surgery which carries potential risks. Please read this leaflet carefully and talk to your doctor about the possible risks.

- Few of the risks are as below**
 - Additional surgery
 - Sometimes, a second surgery may be needed to remove, replace, or adjust the IPL V2.0. It may need to be removed during cataract surgery. Additional surgery might also be needed to reduce eye pressure.
 - Cataract formation
- After Phakic IOL surgery, there is an increased risk that the natural lens of your eye may become cloudy, leading to cataracts that could require surgery. The longer the IPL V2.0 stays in the eye, the higher this risk becomes. Regular eye exams, including IPL V2.0, can cause additional loss of these cells. If too many cells are lost, the cornea can become cloudy and affect your vision. Before surgery, your doctor will examine your eyes to check if you are suitable for IPL V2.0. After the IOL and the cloudy lens are removed, and the long-term is implanted, just like routine cataract surgery.
- Loss of Best Corrected Vision (BCVA)
- The most common reason was clouding of the natural lens (cataract). Vision may get better in these eyes after cataract surgery.
- Raised pressure inside the eye or damage to the optic nerve caused by increased pressure in eye (Glaucoma)
- Normal eye pressure (intraocular pressure or IOP) can vary but is often considered to be from 10-21 millimeters of mercury (mmHg). An IOP higher than normal is called ocular hypertension and if left untreated, can cause damage to the optic nerve and cause permanent vision loss.

Endothelial Cell Loss
The cornea is the clear front part of your eye(s) has a thin layer of special cells called endothelial cells. These cells keep the cornea clear by removing extra fluid. Normally, the number of these cells slowly decreases with age. Eye surgery, including IPL V2.0, can cause additional loss of these cells. If too many cells are lost, the cornea can become cloudy and affect your vision. Before surgery, your doctor will examine your eyes to check if you are suitable for IPL V2.0. After the cornea (called corneal edema), which may require a cornea transplant.

- Other Complications**
 - Other risks associated with IPL V2.0 surgery may include:
 - movement of the colored portion of the eye (iris) through a surgical wound to a position outside the eye (iris prolapse)
 - bleeding in the area on the retina responsible for reading vision (macular hemorrhage),
 - bleeding under the retina (subretinal hemorrhage),
 - increasing in focusing near (accommodation),
 - lifting or pulling of the retina from its natural position (retinal detachment),
 - unequal pupil size (anisocoria).

Adverse reactions and complications due to, or following surgery and implication of any implantable Phakic Contact Lens may include, but are not limited to:

- Increased astigmatism
- Refraction error
- Intraocular infection
- Vitrealous loss
- Visualizing issue
- Anterior subcapsular cataract
- Severe intraocular inflammation
- Glaucoma
- Secondary glaucoma
- Uveitis
- Vitritis
- Hyphema
- Hyphopyon
- Under/over correction
- Non-reactive pupil
- Endophthalmitis
- TASS
- Retinal Detachment
- Iritis
- Cystoid macular edema
- Macular edema

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